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WITH BODY ANALYSIS APP

User Manual

Thanks for choosing the SilverGear Smart Scale. This full body composition scale helps you conveniently measure body weight, body balance, BMI, body fat rate and more.

**IMPORTANT:** Please read the instructions in this manual carefully before using the device. Failure to comply with the instructions or using the scale in ways other than the ones mentioned in this manual may result in injury to yourself or damage to the scale.

#### SPECIFICATION

Product size: 300x260x20mm  
LCD display: 74x28mm  
Weight unit: kg/lb/st  
Weight limit: 6kg-180kg  
Power supply: 2xAAA batteries (not included)  
Account number: each account can add 24 users  
Baby Mode: allows you to scale the weight of a baby  
Height record: allows you to record the height info under 16 years old  
Operating temperature: 10-400C  
Operation humidity: 20%-90%

#### FITDAYS MEASUREMENTS

**Weight:** the weight of the entire body in kg/lb/st.

**Muscle weight:** the weight of the muscles in the body in kg/lb/st.

**Body fat:** the total mass of fat divided by total body mass, multiplied by 100; body fat includes essential body fat and storage body fat.

**Weight without fat:** the weight of the body in kg/lb/st without the weight of body fat.

**Visceral fat:** a type of body fat that wraps around your abdominal organs inside the body.

**Subcutaneous fat:** a type of fat that you develop depending on genetics as well as lifestyle factors such as physical activity and diet.

**Temperature:** the temperature of the body.

**Height record (<16 years old):** the height of children under the age of 16 years.

**Degree of obesity:** an amount of excess weight or body fat that can affect health.

**Calories (daily needs):** calculated daily needs of calorie intake you need in order to get to the weight you entered at your goals.

**Protein (daily needs):** calculated daily needs of protein. Protein is an important building block of bones, muscles, cartilage, skin, and blood.

**Body age:** a measurement of how old you are biologically based upon your health and fitness level as opposed to what your actual age is.

**Bone mass:** the weight of the bones in the body in kg/lb/st.

**Moisture level:** the total amount of fluid that contains in your body. The human body should consist of at least 50% of water.

**BMR:** basal metabolic rate, the minimum number of calories required for basic functions at rest.

**BMI:** body mass index, it combines a person's weight with their height. The results of a BMI measurement can give an idea about whether a person has the correct weight for their height.

**BFR:** blood flow restriction.

#### SUPPORTED OPERATING SYSTEM

iOS 8.0 or higher  
Android 6.0 or higher

#### DOWNLOAD FITDAYS APP

Download "Fitdays" in the App Store or Google Play or scan the QR code:

#### SET UP FITDAYS APP

- Insert 2xAAA batteries at the back of the scale. Place the scale on a hard, flat surface.
- Turn on the Bluetooth and GPS on your smartphone or tablet.
- Register your Fitdays account by email or social media account.
- Fill in your personal information and confirm. Note: Enter correct gender, age and height to get accurate measurements.
- To pair the scale with your smartphone click Account » Device » "+" select the scale ID.

#### IMPORT DATA IN APP

- Tap your foot on the scale to activate the display and shows "0.0"kg.
- Step on the scale to start weighing.  
Note: when stepping on the scale with socks, only weight and BMI will be measured and displayed on the app. To get complete data, step barefoot on the scale.
- The number on the display flashes when the scale is weighing and analyzing.
- When the display shows your left and right weight you have 5-10 seconds to adjust your position, try to get equal left and right weight.
- Continue to stand on the scale until the data is imported in the app. Click Record to see your long-term tracking results.
- Manage users under Account. Add or switch users. The Smart Scale recognizes which user uses the scale and reminds user to switch to the correct user. Each account supports 24 users.
- Baby Mode: click "+" to open Baby Mode. Weigh yourself first then hold the baby to get the weight of the baby.

## TROUBLE SHOOTING GUIDE

### 1. Bluetooth not connected

- Check the operating version of your smartphone, this must be iOS 8.0 or higher/Android 6.0 or higher.
- Click Account, then My Device to check the scale ID.

### 2. No body fat data measured when weighing

- In order to measure body fat stand barefoot on the scale.
- Make sure your feet and the scale are both dry.
- Stand for 10 seconds on the scale. The weighing is finished until the number on the display stops flashing.
- Bluetooth is not enabled.
- Bluetooth is paired with a different scale.

### 3. The scale is producing inconsistent weight

- Check whether the scale is on a hard, flat surface.
- Check each sensor on the back of the scale to make sure nothing sticks to the bottom.
- Remove the batteries to restart the scale.

### 4. The app displays the incorrect fat-muscle ratio

- The body parameter setting is incorrect, check if you entered the correct gender, height, and age.

### 5. How do I re-calibrate the scale after moving it?

- Step on the scale to turn it on. Let it automatically drop to 0.0kg to calibrate.

### 6. The scale does not turn on

- Remove the batteries and put them back. Make sure to insert the batteries in the correct direction.

### 7. How to exchange the unit

- Remove the battery cover, press the unit button to exchange unit from kg to lb or st.

LO	Low battery - please change the batteries (2x1.5V AAA battery)
Err	Overload on the scale (max. 180kg)
C	An error has occurred during measurement. Step off the scale and step onto the scale again to repeat the measuring process.

## WARNING

- Do not use with medical implants such as pacemakers.
- Do not stand on the edge of the scale or jump on it.
- Do not overload the scale (max. 180kg).
- Do not drop or put objects on it as this may damage the sensors.
- Do not immerse the scale in water or use chemical cleaning agents. Clean the scale with slightly damp cloth.
- Not recommended for pregnant women.
- Not recommended for children under 10 years old.
- Any measurement obtained using this device is for reference only and should not be considered as a medical opinion.
- Consult your primary care doctor or physician before making changes to your diet, exercise plan or physical activities.

## BATTERY ADVICE

If the scale isn't used for long periods, it's recommended to remove the battery to avoid damage to the scale due to possible battery leakage.

Keep the battery out of each of children.

Do not mix old and new batteries, with different compositions or of different brands in order to prevent possible leakage.

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